

## BRÜNINGS EXCHANGEABLE TIP

fig. 2; 3.0mm; spoon

for BRÜNINGS interchangeable tubes

**SKU:** 131-486-02

**Categories:** [Goitre](#), [Laryngology](#), [Mouth](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Universal handles and tips</a>
<b>Aesculap</b>	OP196
<b>KLS Martin</b>	39-616-02
<b>AUTHORS</b>	<a href="#">BRUENINGS</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Universal handles and tips</a>
<b>Aesculap</b>	OP196
<b>KLS Martin</b>	39-616-02
<b>AUTHORS</b>	<a href="#">BRUENINGS</a>