

## BRÜNINGS INTERCHANGEABLE TUBE

straight; 22.0cm (8.8")

**SKU:** 131-474-20

**Categories:** [Goitre](#), [Laryngology](#), [Mouth](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Universal handles and tips</a>
<b>Aesculap</b>	OP160
<b>KLS Martin</b>	39-584-20
<b>AUTHORS</b>	<a href="#">BRUENINGS</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Universal handles and tips</a>
<b>Aesculap</b>	OP160
<b>KLS Martin</b>	39-584-20
<b>AUTHORS</b>	<a href="#">BRUENINGS</a>