

BRÜNINGS INTERCHANGEABLE TUBE

straight; 22.0cm (8.8")

SKU: 131-474-20

Categories: [Goitre](#), [Laryngology](#), [Mouth](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|--|
| UNIT | pcs. |
| TYPE | Universal handles and tips |
| Aesculap | OP160 |
| KLS Martin | 39-584-20 |
| AUTHORS | BRUENINGS |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|--|
| UNIT | pcs. |
| TYPE | Universal handles and tips |
| Aesculap | OP160 |
| KLS Martin | 39-584-20 |
| AUTHORS | BRUENINGS |