

## BRUNNER INTESTINAL CLAMP

angled; serrated longitudinal

24.0cm (9.6")

**SKU:** 160-240-24

**Categories:** [Open Surgery](#), [Stomach](#), [Intestines](#), [Rectum](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Stomach and Intestinal Clamp Forceps</a>
<b>Aesculap</b>	EA191
<b>KLS Martin</b>	30-265-24
<b>AUTHORS</b>	<a href="#">BRUNNER</a>
<b>Medicon</b>	50.29.24
<b>SPECIFICATIONS</b>	<a href="#">angled</a> , <a href="#">serrated longitudinal</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Stomach and Intestinal Clamp Forceps</a>
<b>Aesculap</b>	EA191
<b>KLS Martin</b>	30-265-24

**AUTHORS**

[BRUNNER](#)

**Medicon**

50.29.24

**SPECIFICATIONS**

[angled](#), [serrated longitudinal](#)