

BRUNNER LIGATURE CONDUCTOR

strong curved

30.0cm (12")


SKU: 056-500-02

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ligature Conductors
Aesculap	BM932
KLS Martin	20-864-30
AUTHORS	BRUNNER
Medicon	13.18.52
SPECIFICATIONS	curved



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ligature Conductors
Aesculap	BM932
KLS Martin	20-864-30

AUTHORS

[BRUNNER](#)

Medicon

13.18.52

SPECIFICATIONS

[curved](#)