

## BRUNNER LIGATURE CONDUCTOR

strong curved

30.0cm (12")


**SKU:** 056-500-02

**Categories:** [Open Surgery](#), [Suture](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Ligature Conductors</a>
<b>Aesculap</b>	BM932
<b>KLS Martin</b>	20-864-30
<b>AUTHORS</b>	<a href="#">BRUNNER</a>
<b>Medicon</b>	13.18.52
<b>SPECIFICATIONS</b>	<a href="#">curved</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Ligature Conductors</a>
<b>Aesculap</b>	BM932
<b>KLS Martin</b>	20-864-30

**AUTHORS**

[BRUNNER](#)

**Medicon**

13.18.52

**SPECIFICATIONS**

[curved](#)