

BRUNNER LIGATURE CONDUCTOR

left

30.0cm (12")


SKU: 056-500-04

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ligature Conductors
Aesculap	BM934
AUTHORS	BRUNNER
Medicon	13.18.54



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ligature Conductors
Aesculap	BM934
AUTHORS	BRUNNER
Medicon	13.18.54