

BRUNNER LIGATURE CONDUCTOR

slightly curved

30.0cm (12")

SKU: 056-500-01

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ligature Conductors
Aesculap	BM931
KLS Martin	20-863-30
AUTHORS	BRUNNER
Medicon	13.18.51
SPECIFICATIONS	curved



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ligature Conductors
Aesculap	BM931
KLS Martin	20-863-30

AUTHORS

[BRUNNER](#)

Medicon

13.18.51

SPECIFICATIONS

[curved](#)