

BRUNNER LIGATURE NEEDLE

strong curved

24.0cm (9.6")


SKU: 056-455-25

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ligature Needles
Aesculap	BM860
AUTHORS	BRUNNER
SPECIFICATIONS	curved



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ligature Needles
Aesculap	BM860
AUTHORS	BRUNNER
SPECIFICATIONS	curved