

BULLDOG CLAMP ARTERIA; TITAN

strongly curved; 40mm



SKU: 027-045-40

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
MATERIAL	Titanium
KLS Martin	13-145-40
Medicon	15.07.41
SPECIFICATIONS	curved



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
MATERIAL	Titanium
KLS Martin	13-145-40
Medicon	15.07.41
SPECIFICATIONS	curved