

BULLDOG CLAMP ARTERIA; TITAN

strongly curved; 50mm



SKU: 027-045-50

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
MATERIAL	Titanium
KLS Martin	13-145-50
Medicon	15.07.51
SPECIFICATIONS	curved

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
MATERIAL	Titanium
KLS Martin	13-145-50
Medicon	15.07.51
SPECIFICATIONS	curved