

BULLDOG CLAMP ARTERIA; TITAN

straight; 40mm



SKU: 027-042-40

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
MATERIAL	Titanium
KLS Martin	13-142-40
Medicon	15.04.41
SPECIFICATIONS	straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
MATERIAL	Titanium
KLS Martin	13-142-40
Medicon	15.04.41
SPECIFICATIONS	straight