

BULLDOG CLAMP ARTERIA; TITAN

slightly curved; 50mm



SKU: 027-043-50

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bulldog Clamps](#)

MATERIAL [Titanium](#)

KLS Martin 13-143-50

Medicon 15.05.51

SPECIFICATIONS [curved](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bulldog Clamps](#)

MATERIAL [Titanium](#)

KLS Martin 13-143-50

Medicon 15.05.51

SPECIFICATIONS [curved](#)