

BULLDOG CLAMP VENA; TITAN

straight; 40mm



SKU: 027-032-40

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bulldog Clamps](#)

Aesculap FB330

MATERIAL [Titanium](#)

KLS Martin 13-132-40

Medicon 15.04.40

SPECIFICATIONS [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bulldog Clamps](#)

Aesculap FB330

MATERIAL [Titanium](#)

KLS Martin 13-132-40

Medicon 15.04.40

SPECIFICATIONS

[straight](#)