

BULLDOG CLAMP VENA; TITAN

slightly curved; 35mm



SKU: 027-033-35

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bulldog Clamps](#)

MATERIAL [Titanium](#)

KLS Martin 13-133-35

Medicon 15.05.35

SPECIFICATIONS [curved](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bulldog Clamps](#)

MATERIAL [Titanium](#)

KLS Martin 13-133-35

Medicon 15.05.35

SPECIFICATIONS [curved](#)