

BULLDOG CLAMP VENA; TITAN

slightly curved; 40mm



SKU: 027-033-40

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bulldog Clamps](#)

Aesculap FB331

MATERIAL [Titanium](#)

KLS Martin 13-133-40

Medicon 15.05.40

SPECIFICATIONS [curved](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bulldog Clamps](#)

Aesculap FB331

MATERIAL [Titanium](#)

KLS Martin 13-133-40

Medicon 15.05.40

SPECIFICATIONS

[curved](#)