

BULLDOG CLAMP VENA; TITAN

strongly curved; 30mm



SKU: 027-035-30

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
MATERIAL	Titanium
KLS Martin	13-135-30
Medicon	15.07.30
SPECIFICATIONS	curved



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
MATERIAL	Titanium
KLS Martin	13-135-30
Medicon	15.07.30
SPECIFICATIONS	curved