

BUNNELL TENDON STRIPPER

fig. 2; Ø 4.0mm

15.0cm (6")

SKU: 195-268-02

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tendon Stripper
Aesculap	F0932
KLS Martin	23-652-02
AUTHORS	BUNNELL
Medicon	78.20.04

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tendon Stripper
Aesculap	F0932
KLS Martin	23-652-02
AUTHORS	BUNNELL

Medicon

78.20.04