

BUNNELL TENDON STRIPPER

fig. 3; Ø 5.0mm

15.0cm (6")

SKU: 195-268-03

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tendon Stripper
Aesculap	F0933
KLS Martin	23-652-03
AUTHORS	BUNNELL

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tendon Stripper
Aesculap	F0933
KLS Martin	23-652-03
AUTHORS	BUNNELL