

BUNNELL TENDON STRIPPER

fig. 4; Ø 6.0mm

15.0cm (6")

SKU: 195-268-04

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tendon Stripper
Aesculap	FO934
KLS Martin	23-652-04
AUTHORS	BUNNELL
Medicon	78.20.06

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tendon Stripper
Aesculap	FO934
KLS Martin	23-652-04
AUTHORS	BUNNELL

Medicon

78.20.06