

## BUNNELL TENDON STRIPPER

fig. 1; Ø 3.0mm

15.0cm (6")

**SKU:** 195-268-01

**Categories:** [Open Surgery](#), [Orthopaedics](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tendon Stripper</a>
<b>Aesculap</b>	F0931
<b>KLS Martin</b>	23-652-01
<b>AUTHORS</b>	<a href="#">BUNNELL</a>
<b>Medicon</b>	78.20.03

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tendon Stripper</a>
<b>Aesculap</b>	F0931
<b>KLS Martin</b>	23-652-01
<b>AUTHORS</b>	<a href="#">BUNNELL</a>

**Medicon**

78.20.03