

## BUNNELL TENDON STRIPPERS

Ø 5.0mm; straight

23.0cm (9.2")

**SKU:** 195-260-23

**Categories:** [Open Surgery](#), [Orthopaedics](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tendon Stripper</a>
<b>Aesculap</b>	FO940
<b>KLS Martin</b>	23-664-23
<b>AUTHORS</b>	<a href="#">BUNNELL</a>
<b>Medicon</b>	78.26.63
<b>SPECIFICATIONS</b>	<a href="#">straight</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tendon Stripper</a>
<b>Aesculap</b>	FO940
<b>KLS Martin</b>	23-664-23

## AUTHORS

[BUNNELL](#)

## Medicon

78.26.63

## SPECIFICATIONS

[straight](#)