

## BUTTON ELECTRODE; FLEXIBLE

monopolar; 5 Fr.; 470mm

**SKU:** 582-510-05

**Categories:** [Endoscopy](#), [GUBBINI/Gyn](#), [Gynecology](#), [Urology](#)



## ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [HF-Electrodes](#)



## ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [HF-Electrodes](#)