

CAROLL TENDON HOLDING FORCEPS

curved; serrated

11.0cm (4.4")



SKU: 195-245-11

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tendon Seizing Instruments
KLS Martin	23-679-11
AUTHORS	CAROLL
Medicon	78.10.11
SPECIFICATIONS	curved , serrated



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tendon Seizing Instruments
KLS Martin	23-679-11
AUTHORS	CAROLL
Medicon	78.10.11

SPECIFICATIONS

[curved](#), [serrated](#)