

## CAROLL TENDON HOLDING FORCEPS

curved; serrated

11.0cm (4.4")



**SKU:** 195-245-11

**Categories:** [Open Surgery](#), [Orthopaedics](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tendon Seizing Instruments</a>
<b>KLS Martin</b>	23-679-11
<b>AUTHORS</b>	<a href="#">CAROLL</a>
<b>Medicon</b>	78.10.11
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">serrated</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tendon Seizing Instruments</a>
<b>KLS Martin</b>	23-679-11
<b>AUTHORS</b>	<a href="#">CAROLL</a>
<b>Medicon</b>	78.10.11

## SPECIFICATIONS

[curved](#), [serrated](#)