

CARREL HAEMOSTATIC FORCEPS

straight; serrated

8.5cm (3.4")



SKU: 040-010-08

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH100
KLS Martin	13-150-10
AUTHORS	CARREL
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH100
KLS Martin	13-150-10
AUTHORS	CARREL

SPECIFICATIONS

[serrated](#), [straight](#)