

CASPAR DISC RONGEUR

straight; 4 x 14mm; toothed

16.0cm (6.4")



SKU: 240-510-04

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Laminectomy Rongeurs |
| Aesculap | FF534R |
| AUTHORS | CASPAR |
| Medicon | 57.23.34 |
| SPECIFICATIONS | straight , toothed |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------|--------------------------------------|
| UNIT | pcs. |
| TYPE | Laminectomy Rongeurs |
| Aesculap | FF534R |
| AUTHORS | CASPAR |
| Medicon | 57.23.34 |

SPECIFICATIONS

[straight](#), [toothed](#)