

CASTANEDA NEONATAL VASCULAR CLAMP

90° angled; graduated; atraumatic

10.0cm (3.9")




SKU: 142-735-03

Categories: [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vascular Clamps
KLS Martin	24-593-03
AUTHORS	CASTANEDA
SPECIFICATIONS	angled



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vascular Clamps
KLS Martin	24-593-03
AUTHORS	CASTANEDA
SPECIFICATIONS	angled