

## CASTROVIEJO-EPKER CALIPER

straight; effective range 40mm

18.0cm (7.2")



**SKU:** 100-232-40

**Categories:** [Diagnostics](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Measuring instruments</a>
<b>Aesculap</b>	OB792
<b>KLS Martin</b>	35-134-18
<b>AUTHORS</b>	<a href="#">CASTROVIEJO</a> , <a href="#">EPKER</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Measuring instruments</a>
<b>Aesculap</b>	OB792
<b>KLS Martin</b>	35-134-18
<b>AUTHORS</b>	<a href="#">CASTROVIEJO</a> , <a href="#">EPKER</a>