

CASTROVIEJO-EPKER CALIPER

straight; effective range 40mm

9.0cm (3.6")



SKU: 100-231-40

Categories: [Diagnostics](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Measuring instruments
AUTHORS	CASTROVIEJO , EPKER

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Measuring instruments
AUTHORS	CASTROVIEJO , EPKER