

## CATHETER FITTING; CONICAL

Ø 2.5mm; 4mm



**SKU:** 087-919-00

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Stopcocks and Connecting Pieces](#)

**KLS Martin** 18-825-01

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Stopcocks and Connecting Pieces](#)

**KLS Martin** 18-825-01