

## CHERON DRESSING FORCEPS

shanks angled; serrated

25.0cm (10")




**SKU:** 030-303-25

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Sponge and Dressing Forceps</a>
<b>Aesculap</b>	BF092
<b>KLS Martin</b>	14-290-25
<b>AUTHORS</b>	<a href="#">CHERON</a>
<b>Medicon</b>	19.28.75
<b>SPECIFICATIONS</b>	<a href="#">angled</a> , <a href="#">serrated</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Sponge and Dressing Forceps</a>
<b>Aesculap</b>	BF092
<b>KLS Martin</b>	14-290-25

## AUTHORS

[CHERON](#)

## Medicon

19.28.75

## SPECIFICATIONS

[angled](#), [serrated](#)