

## CLAR HEAD LAMP

Ø 55mm; bulb 6 V



**SKU:** 100-150-55

**Categories:** [Diagnostics](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Head Lamps</a>
<b>Aesculap</b>	AC461
<b>KLS Martin</b>	17-517-06
<b>AUTHORS</b>	<a href="#">CLAR</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Head Lamps</a>
<b>Aesculap</b>	AC461
<b>KLS Martin</b>	17-517-06
<b>AUTHORS</b>	<a href="#">CLAR</a>