

CLAR-N LED HEAD LAMP 6V

Ø 55mm; with transformer 230V

SKU: 100-155-00

Categories: [Diagnostics](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Head Lamps
KLS Martin	17-523-04
AUTHORS	CLAR

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Head Lamps
KLS Martin	17-523-04
AUTHORS	CLAR