

CLEANWAVE KERRISON PUNCH; 40° UP CUT

1mm; without ejector

20.0cm (8"); 10mm open; CC-coated




SKU: 240-072-01CC-CW

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Laminectomy Punches |
| MATERIAL | Ceramic(CC) |
| AUTHORS | KERRISON |
| SPECIFICATIONS | 40° up cut , without ejector |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Laminectomy Punches |
| MATERIAL | Ceramic(CC) |
| AUTHORS | KERRISON |
| SPECIFICATIONS | 40° up cut , without ejector |