

## CLEANWAVE KERRISON PUNCH; 40° UP CUT

1mm; without ejector

20.0cm (8"); 10mm open; CC-coated



**SKU:** 240-072-01CC-CW

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

|                       |  |
|-----------------------|--|
| <b>UNIT</b>           | <a href="#">pcs.</a>   |
| <b>TYPE</b>           | <a href="#">Laminectomy Punches</a>                          |
| <b>MATERIAL</b>       | <a href="#">Ceramic(CC)</a>                                  |
| <b>AUTHORS</b>        | <a href="#">KERRISON</a>                                     |
| <b>SPECIFICATIONS</b> | <a href="#">40° up cut</a> , <a href="#">without ejector</a> |



### ZUSÄTZLICHE INFORMATIONEN

|                       |  |
|-----------------------|--|
| <b>UNIT</b>           | <a href="#">pcs.</a>   |
| <b>TYPE</b>           | <a href="#">Laminectomy Punches</a>                          |
| <b>MATERIAL</b>       | <a href="#">Ceramic(CC)</a>                                  |
| <b>AUTHORS</b>        | <a href="#">KERRISON</a>                                     |
| <b>SPECIFICATIONS</b> | <a href="#">40° up cut</a> , <a href="#">without ejector</a> |