

## CLEANWAVE KERRISON PUNCH; 40° UP CUT

1mm; without ejector; thin footplate

18.0cm (7.2"); 10mm open




**SKU:** 240-071-01T-CW

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Laminectomy Punches</a>
<b>MATERIAL</b>	<a href="#">Ceramic(CC)</a>
<b>AUTHORS</b>	<a href="#">CLEANWAVE</a> , <a href="#">KERRISON</a>
<b>SPECIFICATIONS</b>	<a href="#">40° up cut</a> , <a href="#">without ejector</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Laminectomy Punches</a>
<b>MATERIAL</b>	<a href="#">Ceramic(CC)</a>
<b>AUTHORS</b>	<a href="#">CLEANWAVE</a> , <a href="#">KERRISON</a>
<b>SPECIFICATIONS</b>	<a href="#">40° up cut</a> , <a href="#">without ejector</a>