

## CLEANWAVE KERRISON PUNCH; 40° UP CUT

4mm; with ejector; thin footplate

23.0cm (9.2"); 10mm open



**SKU:** 240-073-04AT-CW

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

|                       |   |
|-----------------------|---|
| <b>UNIT</b>           | <a href="#">pcs.</a>                                      |
| <b>TYPE</b>           | <a href="#">Laminectomy Punches</a>                       |
| <b>MATERIAL</b>       | <a href="#">Ceramic(CC)</a>                               |
| <b>AUTHORS</b>        | <a href="#">CLEANWAVE</a> , <a href="#">KERRISON</a>      |
| <b>SPECIFICATIONS</b> | <a href="#">40° up cut</a> , <a href="#">with ejector</a> |



### ZUSÄTZLICHE INFORMATIONEN

|                       |   |
|-----------------------|---|
| <b>UNIT</b>           | <a href="#">pcs.</a>                                      |
| <b>TYPE</b>           | <a href="#">Laminectomy Punches</a>                       |
| <b>MATERIAL</b>       | <a href="#">Ceramic(CC)</a>                               |
| <b>AUTHORS</b>        | <a href="#">CLEANWAVE</a> , <a href="#">KERRISON</a>      |
| <b>SPECIFICATIONS</b> | <a href="#">40° up cut</a> , <a href="#">with ejector</a> |