

CLEANWAVE KERRISON PUNCH; 40° UP CUT

2mm; with ejector; thin footplate

18.0cm (7.2"); 10mm open; CC-coated



SKU: 240-071-02ATCC-CW

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
SPECIFICATIONS	40° up cut, with ejector



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
SPECIFICATIONS	40° up cut, with ejector