

## CLIP APPLIER; ROTATABLE; DETACHABLE

for Ethicon clips medium

Ø 5.0mm; 320mm; for Ø 6.5mm trocars



**SKU:** 612-940-01

**Categories:** [Endoscopy](#), [Laparoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Clip Appliers](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Clip Appliers](#)