

CLIP APPLIER; ROTATABLE; DETACHABLE

for Weck clips small

Ø 5.0mm; 320mm; for Ø 6.5mm trocars



SKU: 612-944-00

Categories: [Endoscopy](#), [Laparoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Clip Appliers](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Clip Appliers](#)