

## CLIP APPLIER; ROTATABLE; DETACHABLE

for Horizon clips small

Ø 5.0mm; 320mm, for Ø 6.5mm trocars



**SKU:** 612-946-00

**Categories:** [Endoscopy](#), [Laparoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Clip Appliers](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Clip Appliers](#)