

## COAKLEY ANTRUM TROCAR

with Luer-Lock tubing connection

Ø 2.3mm; 160mm

**SKU:** 084-190-16

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Trocars</a>
<b>Aesculap</b>	EJ312
<b>KLS Martin</b>	18-515-16
<b>AUTHORS</b>	<a href="#">COAKLEY</a>
<b>Medicon</b>	09.52.33

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Trocars</a>
<b>Aesculap</b>	EJ312
<b>KLS Martin</b>	18-515-16
<b>AUTHORS</b>	<a href="#">COAKLEY</a>

**Medicon**

09.52.33