

COLLIER NEEDLE HOLDER

straight; serrated; fenestrated

12.5cm (5")



SKU: 050-022-12

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM207
KLS Martin	20-332-12
AUTHORS	COLLIER
Medicon	10.04.13
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM207
KLS Martin	20-332-12

AUTHORS

[COLLIER](#)

Medicon

10.04.13

SPECIFICATIONS

[serrated](#), [straight](#)