

## COLLIN HAEMOSTATIC FORCEPS

straight; serrated

14.0cm (5.6")




**SKU:** 040-362-14

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>Aesculap</b>	EA039
<b>KLS Martin</b>	13-268-14
<b>AUTHORS</b>	<a href="#">COLLIN</a>
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>Aesculap</b>	EA039
<b>KLS Martin</b>	13-268-14
<b>AUTHORS</b>	<a href="#">COLLIN</a>

## SPECIFICATIONS

[serrated](#), [straight](#)