

## CONTINUOUS FLOW OPERATING SHEATH

FL:181mm; Ø4.3mm; oval; with 1x stopcock

red; use with: 640-025-22



**SKU:** 640-025-21

**Categories:** [Endoscopy](#), [Gynecology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Hysteroscopy</a>
<b>TYPE</b>	<a href="#">Operating Sheaths</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Hysteroscopy</a>
<b>TYPE</b>	<a href="#">Operating Sheaths</a>