

## CONTINUOUS FLOW OUTER SHEATH

Ø 5.0mm; FL 244mm

white; QUICK-LOCK



**SKU:** 640-015-31

**Categories:** [Endoscopy](#), [Gynecology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Hysteroscopy</a>
<b>TYPE</b>	<a href="#">Diagnostic Sheaths</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Hysteroscopy</a>
<b>TYPE</b>	<a href="#">Diagnostic Sheaths</a>