

CONTINUOUS FLOW OUTER SHEATH

Ø 3.7mm; NL 200mm

QUICK-LOCK



SKU: 640-015-21

Categories: [Endoscopy](#), [Gynecology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Hysteroscopy
TYPE	Diagnostic Sheaths

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Hysteroscopy
TYPE	Diagnostic Sheaths