


COOLEY PEDIATRIC VASCULAR CLAMP

fig. 3; 45° angled; atraumatic

14.0cm (5.5")




SKU: 142-711-03

Categories: [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vascular Clamps
KLS Martin	24-573-03
AUTHORS	COOLEY
SPECIFICATIONS	angled



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vascular Clamps
KLS Martin	24-573-03
AUTHORS	COOLEY
SPECIFICATIONS	angled