

COOLEY SUCTION TUBE

curved; with suction control

inner Ø 6.5mm; outer Ø 8.0mm; WL: 20.0cm

SKU: 085-024-02

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
KLS Martin	18-561-35
AUTHORS	COOLEY
Medicon	09.86.10

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
KLS Martin	18-561-35
AUTHORS	COOLEY
Medicon	09.86.10